



Daily Laboratories

Division of Mobilab, Inc.
2200 W. Altorfer Dr. Suite B · Peoria, IL 61615

ISO 17025
Registry no. 17-B-00133

Using a sponge-on-a-stick procedures

When collecting environmental samples with a sponge on a stick, please remove the entire handle and leave only the sponge in the bag.

While holding the sponge through the bag, the entire handle will twist off. Leaving a portion of the handle attached to the sponge creates sharp edges that causes holes and tears in the outer bag. This creates a cross-contamination opportunity for the samples and a potential biohazard for anyone handling the samples.

Watch here: <https://multimedia.3m.com/mws/media/20162980/3m-fsd-environmental-scrub-sampler-with-stick.mp4>